


Features:

- **Five year parts and one year labor limited warranty***
- Works with Tri-Tech Medical X Series outlets to allow outlets to be installed in walls with 1 3/4" to 2 1/2" of wall covering materials (measured from surface of studs).
- Easy to install – attaches with set screws
- Made in U.S.A.

* See Terms and Conditions, Document No. 99-0477, on our Website at: www.tri-techmedical.com. For complete details.



Ordering Information

Recessed Extension Kit	MRI Compatible
XC1000-EX250	XC1000-EX250-MRI

To place an order or request additional information call toll-free 1-800-253-8692 or fax us at 1-440-937-5060 or e mail us at internationalsales@tri-techmedical.com. Please visit our web-site at www.tri-techmedical.com.