

Submittal Data Sheet




Function:

Male quick-connect fittings are designed to be used on all brands and styles of quick-connect gas outlets. The most common types are shown below. The advantage of the quick-connect male adaptor is the easy snap in/snap out connection it offers.

Features:

- Fast and easy to connect and disconnect
- Color coded, labeled and indexed for specific gas services
- Available in NPT and hose barb styles
- 100 psi maximum working pressure

How to Order:

Adapter	Gas Service	1/8" npt male	1/4" npt male	1/4" hose barb
 Ohmeda™ style male adapter	Oxygen	17-0476	QM-124-4	QM-124-4B
	Medical Air	17-0478	QM-116-4	QM-116-4B
	Medical Vacuum	17-0493	QM-122-4	QM-122-4B
	Nitrous Oxide	QM-104-2	QM-104-4	QM-104-4B
	WAGD	QM-132-2	QM-132-4	QM-132-4B
 Chemetron style male adapter	Oxygen	17-0482	QA-124-4	QA-124-4B
	Medical Air	17-0484	QA-116-4	QA-116-4B
	Medical Vacuum	17-0496	QA-122-4	QA-122-4B
	Nitrous Oxide	QA-104-2	QA-104-4	QA-104-4B
	WAGD	QA-132-2	QA-132-4	QA-132-4B
 Puritan style male adapter	Oxygen	17-0480	QP-124-4	QP-124-4B
	Medical Air	17-0481	QP-116-4	QP-116-4B
	Medical Vacuum	17-0495	QP-122-4	QP-122-4B
	Nitrous Oxide	QP-104-2	QP-104-4	QP-104-4B
	WAGD	QP-132-2	QP-132-4	QP-132-4B