


**Submittal Data Sheet**

**Features and Benefits:**

Med Touch Alarms that have “FUTURE” modules can be expanded or added to with these Expansion Kits. Each kit will include the following:

- **Five-year parts and one-year labor limited warranty\***
- Circuit board programmed for gas ID
- Gas specific transducer to monitor pressures
- Gas specific DISS riser with single check
- DISS x 1/2” tube (to be brazed into line)



TU524TK Kit Shown

\* See Terms and Conditions, Document No. 99-0477, on our website at: [www.tri-techmedical.com](http://www.tri-techmedical.com) for complete details.

**Ordering Information**

Gas Service	Expansion Kit	Transducer	Riser	Gas Board
Nitrous Oxide	TU504TK	DU524	AU8904	<b>35-0029</b> <i>(Programmed gas specific for proper gas ID)</i>
Carbon Dioxide	TU508TK	DU508	AU8908	
Nitrogen	TU512TK	DU512	AU8912	
Medical Air	TU516TK	DU516	AU8916	
Medical Vacuum	TU522TK	DU522	AU8922	
Oxygen	TU524TK	DU524	AU8924A	
WAGD	TU532TK	DU532	AU8932	
Instrument Air	TU52080TK	DU52080	AU892080	

The Gas specific circuit boards are programmed for the gas ID before they are shipped.