


Pressure Transducer, Gauge, and DISS Union Demand Check Assemblies for Medical Gas Manifold Systems

Complete assembly of 4" diameter pressure gauge 0-100 & 0-300 PSI, and gas specific DISS union demand check valve and Adapter fittings. Entire assembly is cleaned for oxygen service.

Note: This assembly is intended to be ordered at the same time as the Tri-Tech Genesys Manifold, to allow retro-fitting of the manifold cabinet at Tri-Tech.



Features:

Unique to Tri-Tech Medical Inc., the PSM assembly

- Tee Assembly using the Tri-Tech Medical Inc. Genesys™ manifold high/low line pressure transducer
- Two DISS union demand check valve
- 4" diameter pressure gauge 0-100 & 0-300 PSI
- Two Brass Adaptors 1/4" x 1/4" NPT & 1/4" x 1/8" NPT
- Brass Plug 1/4" Male NPT
- The transducer which is included with the manifold eliminates the need to purchase an additional high/low source line pressure switch.

How to Order: Easy to use modular ordering system. Fill in 2 blanks to specify the manifold that meets your needs.

PSM - **SSSS** - **GGG**



S = Gas Service

02 = Tri-Gas
12 = Nitrogen
2080 = Instrument Air

G = Gauge Pressure Range

100 psi **300 psi**

with Gas Specific Demand Valve, 1/8" & 1/4" NPT Adapter Fittings for DISS union provided as part of Genesys™ manifold

PSM = Pressure Transducer Manifold Assembly