




**Submittal Data Sheet**

**Project Information**

***Specification***

Y Couplers provide a safe effective means of providing two outlet connections from one source. They include an outlet adapter on one end and two check unit outlets on the other ends.

<b>Y Coupler</b>	<b>Oxygen</b>	<b>Medical Air</b>	<b>Vacuum</b>
 Ohmeda-style, with one Ohmeda-style male adapter x two Ohmeda-style female check unit couplers	YCM-24	YCM-16	YCM-22
 Chemetron-style, with one Chemetron-style male adapter x two Chemetron-style female check unit couplers	YCA-24	YCA-16	YCA-22
 DISS, with one DISS female hand-tight nut and nipple x two DISS male check units	YCD-24	YCD-16	YCD-22