



**Submittal Data Sheet**

**Project Information**

Project Number _____ Approval _____
-------------------------------------

**Specifications**



**Suction Canister  
Ordering Information**

Part Number	Canister Size and Description	Packaging
64-0610	800cc	10/case
64-0612	1200cc	12/case
64-0613	2000cc	12/case
1690	Hydrophobic Filter	each
1955	6' length of 3/16" ID sterile disposable tubing	each
0048250	100' roll of 3/16" ID sterile disposable tubing	100' Roll

Disposable Suction Canisters are designed to be cost-effective, easy and safe to use. Canisters are rigid plastic and completely disposable – simply throw the canister away after each case. Since there is no cleaning, the risk of infection is minimal. Improved set-up, change-out, and suction over liner systems.

Disposable Suction Canisters feature easy-to read critical measurements on frosted panels, also providing a writing area for labeling. Patient port, vacuum port, and pour spout included. Self-sealing lid gives a leak-free seal when suction is applied. Lid assembly includes a hydrophobic filter system housed in a tapered plastic housing. Hydrophobic filter system provides prefilter and final filtration, and protects against overflow. Assortment of canister accessories are available, including plastic and wire rings used to easily mount the suction canisters into slide brackets.

**Slide Bracket Mounting Ring  
Ordering Information**

Part Number	Description	Packaging
64-0620	Plastic ring for 1200cc and 2000cc disposable suction canisters	12/case
64-0621	Wire ring for 1200cc and 2000cc disposable suction canisters	12/case

