



Satin look interchangeable outlet fronts

- Interchangeable with many competitors (Beacon “B”, Amico, Ohio Medical & Gentec) current outlet designs
- Gas specific pin indexing
- Durable scratch resistant vinyl label
- DISS style with fast easy to service Frontall™ front loaded cartridge
- 5 year warranty





How to Order Standard Outlet Fronts Complete ALN91Xx90xx-A

Note: Beacon Nitrous Oxide and WAGD outlets have Beacon specific pin indexing and require the **ALN91Xx90xx-AB** interchangeable fronts

Part Numbers	Gas Service	Connection Style
ALN91XA9004-A	Nitrous Oxide	Allied/ Chemetron
ALN91XA9004-AB	Nitrous Oxide	
ALN91XA9008-A	Carbon Dioxide	
ALN91XA9016-A	Medical Air	
ALN91XA9022-A	Medical Vacuum	
ALN91XA9024-A	Oxygen	
ALN91XA9032-A	WAGD	
ALN91XA9032-AB	WAGD	
ALN91XF9004-A	Nitrous Oxide	DISS
ALN91XF9004-AB	Nitrous Oxide	
ALN91XF9008-A	Carbon Dioxide	
ALN91XF9012-A	Nitrogen	
ALN91XF9016-A	Medical Air	
ALN91XF9022-A	Medical Vacuum	
ALN91XF9024-A	Oxygen	
ALN91XF9032-A	WAGD	
ALN91XF9032-AB	WAGD	
ALN91XF902080-A	Instrument Air	
ALN91XQ9004-A	Nitrous Oxide	Ohmeda/Diamond
ALN91XQ9004-AB	Nitrous Oxide	
ALN91XQ9008-A	Carbon Dioxide	
ALN91XQ9016-A	Medical Air	
ALN91XQ9022-A	Medical Vacuum	
ALN91XQ9024-A	Oxygen	
ALN91XQ9032-A	WAGD	
ALN91XQ9032-AB	WAGD	
ALN91XP9004-A	Nitrous Oxide	Puritan Bennett
ALN91XP9004-AB	Nitrous Oxide	
ALN91XP9016-A	Medical Air	
ALN91XP9022-A	Medical Vacuum	
ALN91XP9024-A	Oxygen	
ALN91XP9032-A	WAGD	
ALN91XP9032-AB	WAGD	
Secondary Check Valve Conversion Kits		
XC9100-C	Beacon B's & Clones	All Positive Pressure Gases
XC9100-CV	Beacon B's & Clones	All Negative Pressure Gases
XC9900CR	Clip Ring Insertion Tool	