



### Tri-Tech look interchangeable outlet fronts

- Interchangeable with many competitors (Beacon & Puritan "B", Amico, Ohio Medical & Gentec) current outlet designs
- Gas specific pin indexing
- Durable scratch resistant vinyl label
- DISS style with fast easy to service Frontall™ front loaded cartridge
- 5 year parts and 1 year labor warranty





**How to Order** Standard Outlet Fronts Complete Xx90xx-A

**Note:** Beacon Nitrous Oxide and WAGD outlets have Beacon specific pin indexing and require the **Xx90xx-AB** interchangeable fronts

Part #	Gas Service	Connection Style
XA9004-A	Nitrous Oxide	Allied/ Chemetron
<b>XA9004-AB</b>	Nitrous Oxide	
XA9008-A	Carbon Dioxide	
XA9016-A	Medical Air	
XA9022-A	Medical Vacuum	
XA9024-A	Oxygen	
XA9032-A	WAGD	
<b>XA9032-AB</b>	WAGD	
XF9004-A	Nitrous Oxide	
<b>XF9004-AB</b>	Nitrous Oxide	
XF9008-A	Carbon Dioxide	
XF9012-A	Nitrogen	
XF9016-A	Medical Air	
XF9022-A	Medical Vacuum	
XF9024-A	Oxygen	
XF9032-A	WAGD	
<b>XF9032-AB</b>	WAGD	
XF902080-A	Instrument Air	
XQ9004-A	Nitrous Oxide	Ohmeda/Diamond
<b>XQ9004-AB</b>	Nitrous Oxide	
XQ9008-A	Carbon Dioxide	
XQ9016-A	Medical Air	
XQ9022-A	Medical Vacuum	
XQ9024-A	Oxygen	
XQ9032-A	WAGD	
<b>XQ9032-AB</b>	WAGD	
XP9004-A	Nitrous Oxide	Puritan Bennett
<b>XP9004-AB</b>	Nitrous Oxide	
XP9016-A	Medical Air	
XP9022-A	Medical Vacuum	
XP9024-A	Oxygen	
XP9032-A	WAGD	
<b>XP9032-AB</b>	WAGD	
<b>Secondary Check Valve Conversion Kits</b>		
XC9100-C	Beacon B's & Clones	All Positive Pressure Gases
XC9100-CV	Beacon B's & Clones	All Negative Pressure Gases
XC9900CR	Clip Ring Insertion Tool	