


Submittal Data Sheet

Features and Benefits:

- **Five year parts and one year labor limited warranty***
- Works with Tri-Tech Medical X Series outlets and slide brackets to allow them to be installed in walls with 1-3/4" to 2-1/2" (without spacer) and 1-1/4" to 1-3/4" (with spacer) of wall covering materials (measured from surface to studs).
- Easy to install – attaches with set screws
- Made in U.S.A.

* See Terms and Conditions, Document No. 99-0477, on our Website at: <https://tri-techmedical.com> for complete details.



Ordering Information

Extension Kits

1-3/4" to 2-1/2" without spacer

Recessed Outlet Extension Kit	MRI Compatible
<i>XC1000-EX250</i>	<i>XC1000-EX250-MRI</i>
Slide Bracket Extension Kit	
<i>XC1144-EX250</i>	

1-1/4" to 1-3/4" with white plastic extension spacer

Recessed Outlet Extension Kit	MRI Compatible
<i>XC1000-EX</i>	<i>XC1000-EX-MRI</i>
Slide Bracket Extension Kit	
<i>XC1144-EX</i>	

White and Pewter aluminum spacers are by special order only and require sufficient quantity. Please contact Tri-Tech Medical to discuss your requirements.