



**Submittal Data Sheet**

**Project Information**



Part Number  
64-0003



64-06CAP-6B-1G



64-06CAP-1/2G



64-05BR

Tri-Tech Medical fluid waste collection bottles are designed with safety and economy in mind. Polycarbonate bottles are durable and reusable. Polycarbonate bottles may be autoclaved or sterilized with sterilizing gases and re-used many times. Polycarbonate bottles permit easy visual inspection of the collection bottle contents.

Collection bottles and cap assemblies are sold separately, and are designed for high air flow capability. The metal cap assemblies are designed with float shut-off, vacuum tube extension, and cap. The metal cap assemblies include a DISS 1220 female wrench tight nut and nipple top connection or a suction tubing nipple connection. The cap assemblies with suction tubing nipple connections feature locking stems and nuts for use with suction bottle and cap assembly mounting brackets. The bottle and cap assembly mounting bracket is made of anodized aluminum, and is precisely designed to mount tightly in slide brackets for a strong even fit.

Part Number	Description
64-0003	1200 cc (1/2 gallon) Polycarbonate bottle only – packaged each or in cases of 6 bottles
64-0004	2800 cc (1 gallon) Polycarbonate bottle only – packaged each or in cases of 6 bottles
64-06CAP-1/2G	Suction cap assembly for ½ gallon bottle with DISS wrench-tight female nut & nipple
64-06CAP-1G	Suction cap assembly for 1 gallon bottle with DISS wrench-tight female nut & nipple
64-06CAP-6B-1/2G	Suction cap assembly for ½ gallon bottle with suction tubing nipple
64-06CAP-6B-1G	Suction cap assembly for 1 gallon bottle with suction tubing nipple
64-05BR	Bottle and cap assembly mounting bracket