



**Submittal Data Sheet**



Part Number - **HT-122-2**



Part Number - **HT-124-2**

***How to Order:***

***Ordering Examples:***

**HT-122-2** = DISS Female Handtight nut & nipple x 1/8" NPT male in Vacuum

<b>Gas Service (Knob Color)</b>	<b>1/8" NPT male</b>	<b>1/4" NPT male</b>	<b>1/4" hose barb</b>	<b>5/16" hose barb</b>
Oxygen (Green)	HT-124-2	HT-124-4	HT-124-4B	
Air (Yellow)	HT-116-2	HT-116-4	HT-116-4B	
Vacuum (White)	HT-122-2	HT-122-4	HT-122-4B	HT-122-5B
Nitrous Oxide (Blue)	HT-104-2	HT-104-4	HT-104-4B	
Nitrogen (Black)	HT-112-2	HT-112-4	HT-112-4B	
WAGD/EVAC (Purple)	HT-132-2	HT-132-4	HT-132-4B	HT-132-5B
Carbon Dioxide (Grey)	HT-108-2	HT-108-4	HT-108-4B	
Instrument Air (Red)	HT-2080-2	HT-2080-4	HT-2080-4B	

DISS female probe adapters connect to gas-specific DISS (Diameter Index Safety System) terminal outlets, and provide a means by which secondary equipment can access gas services. They are made out of chrome plated brass, with color coded hand-tight knobs.

DISS female probe adapters are also available in different configurations and with ISO color code. Please contact Tri-Tech Medical Inc. for those.