

Submittal Data Sheet




Function:

Male quick-connect fittings are designed to be used on all brands and styles of quick-connect gas outlets. The most common types are shown below. The advantage of the quick-connect male adaptor is the easy snap in/snap out connection it offers.

Features:

- Fast and easy to connect and disconnect
- Color coded, labeled and indexed for specific gas services
- Available in NPT and hose barb styles
- 100 psi maximum working pressure

How to Order:

Adapter	Gas Service	1/8" npt male	1/4" npt male	1/4" hose barb
 Ohmeda™ style male adapter	Oxygen	QM-124-2	QM-124-4	QM-124-4B
	Medical Air	QM-116-2	QM-116-4	QM-116-4B
	Medical Vacuum	QM-122-2	QM-122-4	QM-122-4B*
	Nitrous Oxide	QM-104-2	QM-104-4	QM-104-4B
	WAGD	QM-132-2	QM-132-4	QM-132-4B*
 Chemetron style male adapter	Oxygen	QA-124-2	QA-124-4	QA-124-4B
	Medical Air	QA-116-2	QA-116-4	QA-116-4B
	Medical Vacuum	QA-122-2	QA-122-4	QA-122-4B*
	Nitrous Oxide	QA-104-2	QA-104-4	QA-104-4B
	WAGD	QA-132-2	QA-132-4	QA-132-4B*
 Puritan style male adapter	Oxygen	QP-124-2	QP-124-4	QP-124-4B
	Medical Air	QP-116-2	QP-116-4	QP-116-4B
	Medical Vacuum	QP-122-2	QP-122-4	QP-122-4B*
	Nitrous Oxide	QP-104-2	QP-104-4	QP-104-4B
	WAGD	QP-132-2	QP-132-4	QP-132-4B*

* For 5/16" hose barb in Medical Vacuum and WAGD, change the 4B to 5B (example QM-122-5B)

ISO color code and other configurations are available. Please contact Tri-Tech Medical Inc. for those.