



Submittal Data Sheet

Project Information

Tri-Tech Medical fluid waste collection containers are designed with safety and economy in mind. Tough Polycarbonate containers are shatter resistance and durable resistance to implosion. Polycarbonate containers are reusable after cleaning and disinfection or sterilization with sterilizing gases and may be autoclaved many times. Polycarbonate containers permit easy visual inspection of the collection contents with markings.



64-0508 -1500 CC Bottle and 64-0509 2500 CC Bottle

Collection containers comes with easily removable lids with sealing latches and secure leak free vacuum seal, wide opening is designed for ease of cleaning and with built-in slider for use with standard slide brackets eliminates need for separate hanger. The lid assemblies with dependable float design prevents overflow in to tubing and vacuum source when the maximum allowed filling level is reached.



3 different Slide Brackets

*See Terms and Conditions, Document No. 99-0477, on our website at: <https://tri-techmedical.com>. For complete details.

Ordering Information

Part Number	Container Working Volume	Working Temperature	Working Humidity	Working Pressure	Life Span
64-0508	1500 CC	0-100° F	< 80%	0 to -15 psi	3 years
64-0509	2500 CC				
64-0614	Metal Slide Bracket				
64-0615	Die Cast Slide Bracket				
64-0616	Aluminum Slide Bracket				