

## **Female Quick-Connect Couplers**

***Function:***

Female quick-connect coupler fittings are designed for use with all brands and styles of male quick-connect adaptors. The most common types are shown below. The advantage of the quick-connect female coupler is the easy snap in/snap out connection it offers.

***Features:***




- Fast and easy to connect and disconnect
- Color coded, labeled and indexed for specific gas services
- Available in 1/4 M NPT only
- 100 psi maximum working pressure

***How to Order:***

***Ordering Examples:***

QMF-24-4 = Ohmeda style female coupler, Oxygen x 1/4 npt male

QAF-16-4 = Chemetron style female coupler, Air x 1/4 npt male

<b>Adaptor</b>	<b>1/4" npt male</b>
 Ohmeda™ style female coupler	QMF-XX-4
 Chemetron™ style female coupler	QAF-XX-4
 Puritan™ style female coupler	QPF-XX-4

***Substitute "XX" with the following gas services:***

***04 = Nitrous Oxide    16 = Medical Air    22 = Vacuum    24 = Oxygen***